



Status of the POLARIS laser system

Marco Hornung, S. Keppler, H. Liebetrau, A. Kessler, A. Seidel, M. Hellwing, F. Schorcht,
J. Körner, A. Sävert, J. Polz, G. Becker, A.K. Arunachalam, D. Klöpfel, J. Hein, and M. C. Kaluza



Helmholtz Institute Jena, Germany



Institute of Optics and Quantum Electronics, FSU Jena, Germany

- POLARIS overview
- High energy pulse amplification
- High contrast frontend
- Pulse characterization & results
- Conclusion



Petawatt Optical Laser Amplifier for Radiation Intensive Experiments

fully diode-pumped

$\lambda_c = 1030 \text{ nm}$


pulse energy = $6.5 \text{ J} (16.6 \text{ J}) / 4 \text{ J}$ on target

pulse duration = 145 fs

peak intensity = $3.5 \times 10^{20} \text{ W/cm}^2$

temporal contrast for ASE $> 5 \times 10^{12}$

1 shot every 40 s



More than 17000 high energy shots in laser matter experiments !

Used targets:

- thin foils (down to 100 nm thickness)
- water droplets (20 µm)
- sandwich foils
- gas-jet
- hydrogen droplets (cryo-cooled, 6-7 µm) (currently used target)




Diagnostics:


- Multi-channel plate (online, energy, ions)
- Scintillation screen & gateable, fast CCD (online, spatial & energy distribution, ions)
- CR-39 Stack (energy, spatial distribution, ions)
- In-vacuum, soft X-ray camera & bent crystal (spectral & spatial resolution)
- 2 ω or 3 ω online focal spot diagnostics
- target back-side deformation (Nomarski-Interferometer)
- Synchronized side-view (532 nm, 1 ns)
- e⁻ spectrometer
- 150 fs low-energy probe pulse unit,
- variable pre-pulse unit
- backreflection diagnostic for safe laser operation

System Layout

POLARIS is now
Double CPA



Overview: photography's




Control-room for experiments



- POLARIS overview
- High energy pulse amplification
- High contrast frontend
- Pulse characterization & results
- Conclusion


A4 – multipass amplifier

Pump profile:





Output: 8 J @ 0.025

240 J pump energy




Beam profile:



A5 setup

- Currently 120 laser diode stacks $\lambda = 940$ nm, P=2.5kW
- Upgrade to 240 stacks is possible
- Yb³⁺ doped CaF₂ crystal or FP20-Glass with 35 mm p



Size: 6 m x 3 m x 3 m



A5 photography



- Active material:
 - Yb:CaF₂ (1.7 at.%)
 - 65 mm diameter
 - 30 mm thickness
- beam diameter approx. 35 mm (FWHM)
- 9-pass configuration




Yb:CaF₂ made by HellmaMaterials / Jena



A5 pump profile homogenization

Main challenges:

- Homogenize pump profile with 120 Stacks (LabII adaptive strategy ¹)
- Automatisation




Measured pump profile with 780 J pump energy



¹ B. Schmidt, M. Hacker, G. Stobrawa, T. Feurer; LAB2-A virtual femtosecond laser lab, <http://www.lab2.de>



A5 output characteristic

Near-field profile:




16.6 J pulse energy with 10 nm FWHM-bandwidth

Simulation:



Spectral intensity:




- Seed-energy = 2.7 J ($g = 6.1$)
- Energy measured by with *Coherent J-50 MB YAG1535* (verified online with 8%-reflection and another J-50 MB)
- Currently limited by fluence on crystal 2 J/cm^2 (coating-LIDT @ 3.5 J/cm^2)


A. Kessler et al., "16.6 J chirped femtosecond laser pulses from a diode-pumped Yb:CaF₂ amplifier", Optics Letters 2014

A5 low energy (2 Joule) compression & focusing


Pulse duration measurement in the far-field:



Near-field profile:




Far-field profile:



Transport beamline & telescope to compressor:




Spectral intensity:




- POLARIS overview
- High energy pulse amplification
- **High contrast frontend**
- Pulse characterization & results
- Conclusion

Motivation: "old" ASE Characteristics of POLARIS



Measurement of energy, pulse duration and focal spot size.


Double-CPA frontend for nonlinear filtering (XPW)



2 mJ pulse energy in a 130 fs laser pulse.


A1 & A2 – regenerative amplifiers

Beam profile:




A2 output: 50 mJ @ 1 Hz (15 nm FWHM-bandwidth)

Reflective polarisation rotator





- 25 round-trips
- TEM₀₀-mode
- spectral shaping mirrors



Öffner stretcher

- pulses stretched from 90 fs to 20 ps
- hardclip: 45 nm, efficiency 65%



50x50 mm² gold grating, 1200 lines/mm


small footprint: 40x20 cm²



Treacy compressor


- 50x50 mm² gold gratings separated by ≈ 9 cm
- transmission 75%
- small footprint 20x20 cm²

- pulses are compressed to Fourier-limit (Gaussian pulse) of 130 fs
- 2 mJ pulse energy after compression




XPW as nonlinear filter

Experimental setup:





- crossed polarizers with high extinction ratio $R=2 \cdot 10^{-6}$
- 1 m focusing and collimating lenses
- 2 mm BaF₂-crystal
- vacuum-housing due to high intensities




Up to 200 μJ pulse energy!

Near-field profile:



Contrast improvement measurement

- intentionally decreased the ASE contrast by decreasing seed of amplifier !




- POLARIS overview
- High energy pulse amplification
- High contrast frontend
- Pulse characterization & results
- Conclusion

Tiled-Grating Compressor



Stretcher:

- 1480 lines/mm
- 85 fs => 2.6 ns
- GDD = $8.28 \times 10^7 \text{ fs}^2$
- TOD = $-3.71 \times 10^8 \text{ fs}^3$
- FOD = $2.71 \times 10^9 \text{ fs}^4$



Pulse compressor

Tiled-grating:





140 mm beam diameter

M. Hornung et al., Appl Phys B, 101: 93–102, (2010).

Far-field Analysis for tiled-grating alignment

- 2 x specular reflection
- 1 x -1st diffraction order

ideal far field:





Achieved experimental errors:

rotational: < 0.2 μrad

translational: < 20 nm


→ Enough to ensure more than 90 % of the peak intensity for a 140 mm diameter, f/2-focussed, 150-fs-pulse (compared to a monolithic grating).

Piston π :




Tiled-grating alignment & shortest pulses

Spectrally resolved far-field measurements:




M. Hornung *et al.*, Opt. Lett., **35** 2073, (2010).

SRFF measurement setup:




Pulse duration measurement:





High-dynamic pulse duration measurements

Installation of a Wizzler & Dazzler loop:




A4 pulse duration:




Problem: 10 year old roof-mirror



Installation of a Wizzler & Dazzler loop



Degraded 10 year old roof-mirror:




Improved ps-contrast:



Further contrast improvement: SHG of the laser pulses!

Experiment aimed in autumn with 160 mm diameter KDP-crystal.

$$I_{\text{Peak}} = 3.5 \times 10^{20} \text{ W/cm}^2$$



ASE at 2×10^{-13} relative intensity (38nJ) with XPW frontend.

Picosecond Contrast Measurements



For 150 fs-Pulses:

2×10^{-13} @ -200 ps (ASE)

1×10^{-9} @ -25 ps


1×10^{-6} @ -5 ps

Old measurement:



To do:
Some post/pre-pulses (e.g. @ -130 ps, 10^{-7}) has to be eliminated in the future!




TNSA p⁺-acceleration: target thickness scan



Target thickness for highest p⁺-energy decreases with improved contrast!


Experimental Performance - TNSA

Maximum proton cut-off energy of 17.7 MeV with 2.3 J of laser energy and a 2 μm Ta-foil.



Laser performance



Focal Spot in Targetarea with adaptive optics:



$t_{FWHM} = 145 \text{ fs}$
 $q = 0.46$
 $A_{FWHM} = 7.1 \mu\text{m}^2$
 $E_{pulse} = 6.5 \text{ Joule}$
 $\eta_{Compressor} = 0.63$

$$I_{\text{Peak}} = 3.5 \times 10^{20} \text{ W/cm}^2$$

Pulse duration:



POLARIS as a fully diode-pumped laser is in experimental operation.

(more than **17000** high intensity experiment-shots)

Improvements in terms of temporal contrast and pulse energy was done.

Pulse parameters:

- 3.5×10^{20} W/cm² peak intensity
- 4 J on target pulse energy (16 J out of A5 soon available)
- 145 fs pulse duration
- Relative temporal contrast of 2×10^{-13} for ASE and 10^{-8} for (some residual) pre-pulses.

M. Hornung et al., "High intensity, high contrast laser pulses generated from the fully diode-pumped Yb:glass laser system POLARIS", Optics Letters 2013.
A. Kessler et al., "16.6 J chirped femtosecond laser pulses from a diode-pumped Yb:CaF₂ amplifier", Optics Letters 2014

Thank You!

